Customer Pro-Strain Assembly Bill of Materials

	٠.	
ĸ	ıt	
ı	ı	.Э

	<u>, </u>			Kits		
		CXPS240001	CXPS260001	CXPS320001	CXST240001	CXST300001
		Pro-Strain	Pro-Strain	Pro-Strain	Pro-Strain	Pro-Strain
		straining/turning	straining/turning	straining/turning	strut kit 2400	strut kit
		post kit 2400 x	post kit 2600 x	post kit 3200 x	x 48 x 3.2mm	3000 x 48 x
Posts	Muti-Packs	89 x 3mm	89 x 3mm	89 x 3mm c/w		3.2mm
	Pro-Strain straining/turning post 2400 x 89 x					
CXPS240002	3mm (24 units) **C1**	1				
	Pro-Strain straining/turning post 2600 x 89 x					
CXPS260002	3mm (24 units) **C2**		1			
	Pro-Strain straining/turning post 3200 x 89 x					
CXPS320002	3mm (24 units) **C3**			1		
CXPS100001	Pro-Strain strut peg 1000 x 89 x 3mm (24 units)				1	1
Struts						
CXST240002	Pro-Strain strut 2400 x 48 x 3.2mm (24 units)				1	
CXST300002	Pro-Strain strut 3000 x 48 x 3.2mm (24 units)					1
Plates*						
	Pro-Strain breast/heel plate 250 x 250 x 4mm					
CXFX250001	(10 units)	2	2	2	1	1
	Pro-Strain breast/heel plate 350 x 250 x 4mm					
CXFX350001	(10 units)	2*	2*	2*	1*	1*
	Pro-Strain breast/heel plate 600 x 250 x 4mm (4					
CXFX600001	units)	2#	2#	2#	2#	2#
Strut adjuster			,	,		
CXFX000001	Pro-Strain adjuster and coupler set (10 units)				1	1
Post cap						
CXFX000002	Pro-Strain post caps 89mm (10 units)	1	1	1	1	1
Nuts & bolts						
CXFX000003	Pro-Strain fixings M12 x 110mm (100 units)	4	4	4	2	2

^{*} Large plate to be used in soft ground

Note: Box assembly replace the 'strut peg' with another straining post. A Gripple G-pak for bracing wire will also be required

Extra-large plate often used on strut pets in conjunction with a 250mm plate as per the illustration

Note: Half box is required for Deer fencing use a 2.4m strainer and the secondary post.